



CO/6100

Decon 90



The leading surface active cleaning agent and radioactive decontaminant, for laboratory, medical and industrial critical cleaning applications. Biodegradable, bactericidal, phosphate free, non-flammable and totally rinsable. Offering improved performance without the unacceptable hazards associated with the storage, handling, use and disposal of corrosive acids. Decon 90 provides consistent and reliable cleaning and decontamination in the modern environment.

Suitable for cleaning glassware, plastics, ceramic surfaces and ferrous metals.

Cat. No.	Size, litre
CD/61002	1
CD/61004	5
CD/61008	20



CO/6130

Neutracon



A specialised surface active cleaning agent for safe cleaning and decontamination of sensitive materials which may be corroded, etched, discoloured or weakened by acidic or alkaline cleaning agents. Biodegradable, phosphate free, non-corrosive, non-flammable and totally rinsable. Specifically formulated to provide an effective near neutral cleaning agent with maximum cleaning and decontaminating properties.

Decon Neutracon is safe to use on non-ferrous metals, soft glass, coated lenses, polymer surfaces and other sensitive materials which must not be cleaned by acid or alkaline cleaning agents.

Cat. No.	Size, litre
CD/61302	1
CD/61304	5



CO/6140

Dri-decon



A highly effective surface active cleaning agent, in granular form, for use in all types of glassware washing machines. Biodegradable, non-flammable, non-foaming and totally rinsable. Dri-decon is a highly effective cleaning and decontaminating agent specially formulated to give maximum performance with all types of glass washing machines, which use a main wash detergent in 'powder' form.

Cat. No.	Size, kilo
CD/61404	5



CO/6110

Decon Decomatic



A liquid surface active cleaning agent for use in all types of laboratory, medical and industrial glassware washing machines. Biodegradable, bactericidal, phosphate free, non-foaming, non-flammable and totally rinsable. Decon Decomatic has been formulated to provide outstanding detergent and free rinsing action for improved results from every type of glassware washing machine which accepts a main wash detergent in liquid form. May also be used as a pre-soak prior to machine washing.

Cat. No.	Size, litre
CD/61104	5
CD/61106	10



CO/6120

Decon Acid Rinse



Acid rinse (neutraliser) and pre-wash (preliminary) cleaner. For use in all types of laboratory, medical and specialised industrial glassware washing machines used in conjunction with main wash detergents. Biodegradable, bactericidal, non-foaming, non-flammable and totally rinsable. Decon Acid Rinse improves the performance of glassware washing machines by preventing the carry-over of traces of the main wash detergent which may be alkaline or contain phosphates and by inhibiting the deposition of salts from 'hard' water supplies.

Cat. No.	Size, litre
CD/61204	5

For Neodisher glassware washing machine detergents, please see page 515.