



MS/16202

Microscope slide staining jar

Robust polypropylene jars, can hold up to 10 76 x 26mm slides back to back. Can also be used as a mailing container. Supplied complete with lid.

Cat. No.	Code	Description
MS/16202	SWM012	Staining jar, square pattern



MS/92116

Microscope slide staining trough

Glass trough with loose flat lid 85 x 105 x 75mm high. Accepts one glass accessory tray which has inner grooves to accept 10 slides 76 x 26mm horizontally back to back. Must be ordered separately.

Cat. No.	Description
MS/92116	Staining trough
MS/92118	Accessory tray



MS/92204

Microscope slide staining trough simple, Hellendahl

Glass trough with loose flat lid. Internally grooved to accept 16 slides 76 x 26mm vertically back to back

Cat. No.	Description
MS/92204	Staining trough



MS/92206

Microscope slide staining trough simple, Schiefferdecker

Glass trough with loose flanged lid. Internally grooved to accept 10 slides 76 x 26mm horizontally back to back

Cat. No.	Description
MS/92206	Staining trough



MS/92208

Microscope slide staining trough, square, large

Glass trough with drop-on-lid. Approximate int. dims 115 x 115 x 77mm.

Cat. No.	Description
MS/92208	Staining trough