



CV/6327

## Chromatography vials (continued)

Cat. No.	Code	Description	Pack Qty
CV/63240	5400SB-08	8mm Clear LDPE Snap Plug	200
CV/63242	5400SB-08R	8mm Red LDPE Snap Plug	200
CV/63244	5400SB-08B	8mm Blue LDPE Snap Plug	200
CV/63246	5400SB-08G	8mm Green LDPE Snap Plug	200
CV/63248	5400SB-08Y	8mm Yellow LDPE Snap Plug	200
CV/63250	5405SB-12	12mm PE Starburst Conical Snap Plug Opaque/natural PE *	100
CV/63252	601010-08	8mm White PTFE Septa	1000
CV/63254	601010-13	13mm White PTFE Septa	1000
CV/63256	601080-08	White PTFE with slit, 0.25mm thick	1000
CV/63258	603090-08	Black Viton, 0.7mm thick	1000
CV/63260	604040-08	PTFE/Butyl rubber, 1.20mm thick	1000
CV/63262	604060-08	Red PTFE/Silicone, 1.00mm thick	1000
CV/63264	606040-13	PTFE/Rubber 1.2mm	1000
CV/63266	606070-08	Red PTFE/Silicone with slit 1.50mm thick	1000
CV/63268	604075-08	Red PTFE/Silicone/Red PTFE with slit, 1.00mm thick	1000
CV/63270	604550-08	Red PTFE/Silicone with slit 1.20mm thick	1000
CV/63272	606550-08	Red PTFE/Silicone 1.65mm thick	1000
CV/63274	606550B-08	Blue PTFE/Silicon 1.65mm thick	1000
CV/63276	606550T-08	Tan PTFE/Silicone 1.65mm thick	1000
CV/63278	607550-08	Red PTFE/Silicone, 1.90mm thick	1000
CV/63280	610050-08	PTFE/Silicone, 2.50mm thick	1000

\* Other colours available, please ask for details.



CD/12000

## Cleaning

## Lipsol Detergent

- Excellent cleaning performance, ideal for use with laboratory glassware, plasticware, rubber, ceramics and stainless steel
- Economical price
- Non-corrosive, non-irritant and totally biodegradable surface-active cleaning agent
- Free of chlorine bleaches, phosphates and enzymes

Cat No.	Description	Pack Qty
CD/12000	Lipsol concentrate	2 x 5 Litre Jerry Cans

## Disinfectant, virucidal, Virkon

A wide spectrum disinfectant which has been proven to be effective against all 22 virus families affecting man and animals, including HIV and Hepatitis B, and successfully tested on over 400 strains of bacteria. In short, Virkon kills all known key human pathogens.

Virkon is still the ultimate high level surface disinfectant suitable for use on surfaces in Healthcare, Laboratory, Food preparation and production, Catering, Residential, Hotel, Transport, Emergency Services and Public environments. In addition to Virkon powder a new 5gm tablet provides a quick way to make up 500ml of 1% Virkon solution. The tablet dissolves cleanly and quickly, and is conveniently shaped to pop into the neck of a trigger spray bottle for rapid and safe application of the solution.

Virkon can also be applied, in powder form, direct to spillages of blood and other body fluids. Non-corrosive when used as directed, non-bleaching, bio-degradable. Use 1% solution.



CO/5100

Cat. No.	Size	Pk qty
CD/51002	50g sachet	50
CD/51004	500g squeeze "puff" pack	6
CD/51006	1 x 5kg drum	each
CD/51008	10 x 5 gm Tablets	20 tubs
CD/51010	50 x 5 gm Tablets	12 tubs